

## Product datasheet for **RC221780L2V**

### TROVE2 (RO60) (NM\_004600) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TROVE2 (RO60) (NM_004600) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RO60
Synonyms:	RORN; SSA2; TROVE2
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004600
ORF Size:	1614 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221780).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004600.3</a>
RefSeq Size:	3098 bp
RefSeq ORF:	1617 bp
Locus ID:	6738
UniProt ID:	<a href="#">P10155</a>
Cytogenetics:	1q31.2
Protein Families:	Transcription Factors
Protein Pathways:	Systemic lupus erythematosus



[View online »](#)

**MW:** 60.5 kDa

**Gene Summary:** RNA-binding protein that binds to misfolded non-coding RNAs, pre-5S rRNA, and several small cytoplasmic RNA molecules known as Y RNAs. May stabilize some of these RNAs and protect them from degradation (PubMed:18056422). Binds to endogenous Alu retroelements which are induced by type I interferon and stimulate proinflammatory cytokine secretion. Regulates the expression of Alu retroelements as well as inflammatory genes (PubMed:26382853).[UniProtKB/Swiss-Prot Function]