

Product datasheet for **RC221805A1V**

Human C3orf32 (SSUH2) (NM_015931) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human C3orf32 (SSUH2) (NM_015931) AAV Particle
Tag:	Myc-DDK
Symbol:	C3orf32
Synonyms:	C3orf32; fls485; SSU-2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC221805 representing NM_015931 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCCCTCCCCAGTTGGACTGCTGAGGGCTTTACCACTACCGTGGCCTCAGTTTCTCGCCTGCACGTTGA
 GGAGGCTGGCTGGCCCGCTGAGTCCACAGGCCCTTCCAGAAGCCCCGCTCTCTGTTCGGTCCCCTG
 CAGAGTCCCTGCGATGACGGAGGAGTGGCCCGGAAGCCCTCCTCAGCTTTGTGGACTCTAAGTGCTGC
 TACAGCAGCACGGTGGCTGGAGACCTCGTCATCCAGGAGCTGAAGCGGCAGACCCTCTGCAGGTACCGTC
 TGGAGACCTTTAGTGAATCCAGGATAAGCGAGTGGACATTTCAACCCTTTACTAACCACCTCTGTGGATGG
 GCCGCAAGAGGCGCCTCCCCAGGCTCTGGGACATCAAGGTTCAAGGTCCTCCGATGTTTCAGGAAGAC
 ACCAGGAAGTTCAGGTCCCTCACTCGTCACTGGTCAAGGAATGCCACAAATGCCATGGGCGTGGGCGGT
 ACAAGTGCAGCGGCTGCCACGGGGCGGCACGGTGCGGTGCCATCCTGCTGCGGAGCCAAGCGCAAAGC
 CAAGCAGTCCCGGAGATGTCAGCTGTGCGCGGGTCCGGCAGGCGAAGATGCAGCACTTGCTCAGGGAGA
 GGGAACAAGACCTGCCACCTGCAAGGGGAGAAGAAGCTGTTGCACTTCATCCAGCTTGTCATCATGT
 GGAAGAACAGCTTGTGAGTTGTGTCTGAGCACCGGCTCAACTGCCCCAGGGAGCTCCTTGCTAAAGC
 CAAAGGAGAAAACCTCTTAAGGATGAAAACCTCGGTGGTGTACCCCATCGTGGACTTCCCTCTGCGAGAC
 ATCTCTCTTGCTCCAGAGGGGCATTGCAGAGCACAGCGCTGCCTTGGCCTCCCGTGCCCGCTCTGCG
 AGCAGCGCCAGACCATGAGCTGATCCCTCACAAGGTTCACTATTGGTACCAAGGAAAGACTTATGT
 CTACTACATCTATGGCACTGACCACCAGGTGTATGCGGTGGACTATCCTGAGCGGTATTGCTGTGGCTGT
 ACCATCGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>RC221805 representing NM_015931 Red=Cloning site Green=Tags(s) MPSPVGLLRALPLPWPQFLACTLRRLAGPRESTGPSQKPPPLCSVPCRPAMTEEVAREALLSFVDSKCC YSSTVAGDLVIQELKRQTLCRYRLETFSERISEWTFQPFTNHSVDGPQRGASPLWDIKVQGPPMFQED TRKFQVPHSSLVKECHKCHGRGRYKCSGCHGAGTVRCPSCCGAKRKAKQSRRQCAGSGRRRCSTCSGR GNKTCATCKGEKKLLHFIQLVIMWKNSLFEFVSEHRLNCPRELLAKAKGENLFKDENSVVYPIVDFPLRD ISLASQRGIAEHSAAASRARVLQQRQTIELIPLTEVHYWYQGKTYVYIYGTDHQVYAVDYPERYCCGC TIV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_015931
ORF Size:	1059 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_015931.1
RefSeq Size:	1545 bp
RefSeq ORF:	1062 bp
Locus ID:	51066
UniProt ID:	Q9Y2M2
Cytogenetics:	3p25.3
MW:	39.7 kDa