

Product datasheet for **RC221856L3V**

DLGAP1 (NM_001003809) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DLGAP1 (NM_001003809) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DLGAP1
Synonyms:	DAP-1; DAP-1-ALPHA; DAP-1-BETA; DAP1; DLGAP1A; DLGAP1B; GKAP; SAPAP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001003809
ORF Size:	2025 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221856).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001003809.1 , NP_001003809.1
RefSeq Size:	2393 bp
RefSeq ORF:	2028 bp
Locus ID:	9229
UniProt ID:	O14490
Cytogenetics:	18p11.31
MW:	75.1 kDa
Gene Summary:	Part of the postsynaptic scaffold in neuronal cells.[UniProtKB/Swiss-Prot Function]



[View online »](#)