

## Product datasheet for **RC221868L3V**

### CLEC4F (NM\_173535) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CLEC4F (NM_173535) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CLEC4F
Synonyms:	CLECSF13; KCLR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_173535
ORF Size:	1767 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221868).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_173535.2</a>
RefSeq Size:	2486 bp
RefSeq ORF:	1770 bp
Locus ID:	165530
UniProt ID:	<a href="#">Q8N1N0</a>
Cytogenetics:	2p13.3
Protein Families:	Druggable Genome, Transmembrane
MW:	65.3 kDa



[View online »](#)

**Gene Summary:**

Receptor with an affinity for galactose and fucose. Could be involved in endocytosis (By similarity).[UniProtKB/Swiss-Prot Function]