

Product datasheet for **RC221869**

PRELID2 (NM_205846) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PRELID2 (NM_205846) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: PRELID2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC221869 representing NM_205846
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGGGGTCTCGGTGGATGTGCACCAGGTGTACAAGTACCCCTTCGAGCAGGTGGTCCGAGCTTTCTCC
GAAAGTACCCCAACCCCATGGATAAAAAATGTCATCTCAGTAAAAATCATGGAGGAAAAAGAGATGAATC
AACAGGGGTCACTACAGAAAGAGGATTGCAATCTGTCAGAACGTGGTCCAGAAATTTAAGGAAGGTG
AGCATTTTGAAGTACCTAATATCCAATTAGAAGAGGAGTCATGGCTCAATCCTCGGAAAGAAACATGG
CCATACGGAGTCACTGCCTTACGTGGACACAGTATGCATCCATGAAGGAAGAGTCTGTCTCCGGGAAAG
TATGGAAAACCAAATTGGACAGAGTTCATTCAAAGAGGCAGGATTTCAATCACAGGGGTTGGATTTCTC
AACTGTGTTTTAGAACTTTTGGCAGCACATTCTACGACAGGGAGCCAGAAAGGAATTAGAATCATGG
AGATGCTGCTAAAGGAACAGTGTGGTGCCTTCTAGCTGAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC221869 representing NM_205846
Red=Cloning site Green=Tags(s)

MGVSDVHQVYKYPFEQVVASFLRKYPNPMKDNVISVKIMEEKRDESTGVIYRKRIAICQNVVPEILRKV
SILKVPNIQLEESWLNPRERNMAIRSHCLTWTQYASMKEESVFRESMENPNWTEFIQRGRISITGVGFL
NCVLETFASTFLRQGAQKIRIMEMLLKEQCGAPLAE

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_205846

ORF Size: 531 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

| | |
|-------------------------------|---|
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_205846.3</u> |
| RefSeq Size: | 2104 bp |
| RefSeq ORF: | 534 bp |
| Locus ID: | 153768 |
| UniProt ID: | <u>Q8N945</u> |
| Cytogenetics: | 5q32 |
| MW: | 20.4 kDa |