

Product datasheet for **RC221884L2V**

CK1 epsilon (CSNK1E) (NM_001894) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | CK1 epsilon (CSNK1E) (NM_001894) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | CK1 epsilon |
| Synonyms: | CK1ε; CK1εpsilon; HCK1E |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_001894 |
| ORF Size: | 1248 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC221884). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001894.4 , NP_001885.1 |
| RefSeq Size: | 2670 bp |
| RefSeq ORF: | 1251 bp |
| Locus ID: | 1454 |
| UniProt ID: | P49674 |
| Cytogenetics: | 22q13.1 |
| Domains: | pkinese, TyrKc, S_TKc |
| Protein Families: | Druggable Genome, Protein Kinase |



[View online »](#)

Protein Pathways: Circadian rhythm - mammal, Hedgehog signaling pathway, Wnt signaling pathway

MW: 47.3 kDa

Gene Summary: The protein encoded by this gene is a serine/threonine protein kinase and a member of the casein kinase I protein family, whose members have been implicated in the control of cytoplasmic and nuclear processes, including DNA replication and repair. The encoded protein is found in the cytoplasm as a monomer and can phosphorylate a variety of proteins, including itself. This protein has been shown to phosphorylate period, a circadian rhythm protein. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Feb 2014]