

Product datasheet for **RC221918L3V**

EHBP1L1 (NM_001099409) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EHBP1L1 (NM_001099409) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EHBP1L1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001099409
ORF Size:	4569 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221918).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001099409.1
RefSeq Size:	5197 bp
RefSeq ORF:	4572 bp
Locus ID:	254102
UniProt ID:	Q8N3D4
Cytogenetics:	11q13.1
MW:	161.7 kDa
Gene Summary:	May act as Rab effector protein and play a role in vesicle trafficking.[UniProtKB/Swiss-Prot Function]



[View online »](#)