

Product datasheet for **RC221932L4V**

SLC8A3 (NM_058240) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC8A3 (NM_058240) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC8A3
Synonyms:	NCX3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_058240
ORF Size:	2772 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221932).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_058240.2
RefSeq Size:	5259 bp
RefSeq ORF:	2775 bp
Locus ID:	6547
UniProt ID:	P57103
Cytogenetics:	14q24.2
Domains:	Calx-beta, Na_Ca_Ex
Protein Families:	Transmembrane



[View online »](#)

Protein Pathways: Calcium signaling pathway

MW: 99.1 kDa

Gene Summary: This gene encodes a member of the sodium/calcium exchanger integral membrane protein family. Na⁺/Ca²⁺ exchange proteins are involved in maintaining Ca²⁺ homeostasis in a wide variety of cell types. The protein is regulated by intracellular calcium ions and is found in both the plasma membrane and intracellular organellar membranes, where exchange of Na⁺ for Ca²⁺ occurs in an electrogenic manner. Alternative splicing has been observed for this gene and multiple variants have been described. [provided by RefSeq, Aug 2013]