

Product datasheet for **RC221939L1V**

C1GALT1 (NM_020156) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C1GALT1 (NM_020156) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C1GALT1
Synonyms:	C1GALT; T-synthase
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_020156
ORF Size:	1089 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221939).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020156.1
RefSeq Size:	1794 bp
RefSeq ORF:	1092 bp
Locus ID:	56913
UniProt ID:	Q9NS00
Cytogenetics:	7p22.1-p21.3
Protein Families:	Transmembrane
Protein Pathways:	Metabolic pathways, O-Glycan biosynthesis



[View online »](#)

MW: 42 kDa

Gene Summary: The protein encoded by this gene generates the common core 1 O-glycan structure, Gal-beta-1-3GalNAc-R, by the transfer of Gal from UDP-Gal to GalNAc-alpha-1-R. Core 1 is a precursor for many extended mucin-type O-glycans on cell surface and secreted glycoproteins. Studies in mice suggest that this gene plays a key role in thrombopoiesis and kidney homeostasis. [provided by RefSeq, Sep 2010]