

## Product datasheet for **RC221947L3V**

### CPLX2 (NM\_001008220) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CPLX2 (NM_001008220) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CPLX2
Synonyms:	921-L; CPX-2; CPX2; Hfb1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001008220
ORF Size:	402 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221947).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001008220.1</a>
RefSeq Size:	4676 bp
RefSeq ORF:	405 bp
Locus ID:	10814
UniProt ID:	<a href="#">Q6PUV4</a>
Cytogenetics:	5q35.2
Protein Families:	Druggable Genome
MW:	15.4 kDa


[View online »](#)

**Gene Summary:**

Proteins encoded by the complexin/synaphin gene family are cytosolic proteins that function in synaptic vesicle exocytosis. These proteins bind syntaxin, part of the SNAP receptor. The protein product of this gene binds to the SNAP receptor complex and disrupts it, allowing transmitter release. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]