

Product datasheet for RC221955L2

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OVCA1 (DPH1) (NM_001383) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: OVCA1 (DPH1) (NM_001383) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: OVCA1

Synonyms: DEDSSH; DPH2L; DPH2L1; OVCA1

Mammalian Cell None

Selection:

Vector: pLenti-C-mGFP (PS100071)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC221955).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_001383

ORF Size: 1329 bp





OVCA1 (DPH1) (NM_001383) Human Tagged Lenti ORF Clone - RC221955L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001383.3</u>

RefSeq Size: 2200 bp
RefSeq ORF: 1317 bp
Locus ID: 1801

UniProt ID: Q9BZG8
Cytogenetics: 17p13.3

Domains: Diphthamide syn

MW: 48.6 kDa


Gene Summary: The protein encoded by this gene is an enzyme involved in the biosynthesis of diphthamide, a

modified histidine found only in elongation factor-2 (EEF2). Diphthamide residues in EEF2 are targeted for ADP-ribosylation by diphtheria toxin and Pseudomonas exotoxin A. Defects in this gene have been associated with both ovarian cancer and autosomal recessive intellectual disability with short stature, craniofacial, and ectodermal anomalies. [provided by RefSeq, Oct

2016]



Product images:



Circular map for RC221955L2