

Product datasheet for **RC221999L4V**

SEL1L (NM_005065) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SEL1L (NM_005065) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SEL1L
Synonyms:	Hrd3; PRO1063; SEL1-LIKE; SEL1L1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005065
ORF Size:	2382 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC221999).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005065.3
RefSeq Size:	7885 bp
RefSeq ORF:	2385 bp
Locus ID:	6400
UniProt ID:	Q9UBV2
Cytogenetics:	14q31.1
Protein Families:	Druggable Genome, Transmembrane
MW:	88.6 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is part of a protein complex required for the retrotranslocation or dislocation of misfolded proteins from the endoplasmic reticulum lumen to the cytosol, where they are degraded by the proteasome in a ubiquitin-dependent manner. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]