

Product datasheet for **RC222000**

MARCH I (MARCH1) (NM_017923) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: MARCH I (MARCH1) (NM_017923) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: MARCH I
Synonyms: MARCH-I; MARCH1; RNF171
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC222000 representing NM_017923
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACCAGCAGCCACGTTTGTGTAATTTTTAAATATGTGGAAAAAAGCAAGATATCAACCATGTATT
ACCTTAACCAAGATGCCAAATTCTAACTTGTCTCCAGGCAAGCAGCCCAACAACAGGGACAGCTCC
CAGGAGCCAGTCAAGTTGTCTGTCTGTCCATCCACTCAGGACATCTGCAGAATCTGCTACTGCGAAGG
GATGAAGAGAGCCCCCTCATCACACCCTGTCGCTGCCTGGGACTGCGCTTTGTCCACCAGTCTGCC
TCCACCAGTGGATAAAGAGCTCAGATACACGCTGCTGTGAGCTCTGCAAGTATGACTTCATAATGGAGAC
CAAGCTCAAACCCCTCCGGAAGTGGGAGAACTACAGATGACCACAAGTAAAGGAGGAAAAATTTCTGC
TCTGTACATTCCACGTAATCGCGATCACCTGTGTGGTTTGGTCTTTGTATGTATTGATAGACCGGACAG
CGGAGGAAATCAAGCAAGGCAATGACAATGGTGTCTTGAATGGCCATTTTGGACAAAACTGGTTGTGGT
AGCCATTGGCTTACAGGAGGTCTTGTCTTACGTACGTACAGTAAAGTCTATGTTTCAGTTGTGGCGC
AGGCTGAAGGCCACAACCGTGTGATCTTTGTACAAAATTGCCAGACTGCCAAAACTGGAGAAGA
ACTTCTCATGTAATGTAACACAGACATCAAAGATGCTGTGGTAGTGCCTGTACCACAAACAGGTGCAAA
TTCCTGCCATCTGCAGAGGGTGGCCCCCTGAAGTTGTATCAGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC222000 representing NM_017923
 Red=Cloning site Green=Tags(s)


MTSSHVCCNFLNMWKKSKISTMYLNDQAKLSNFLQASSPTTGTAPRSQSRLSVC PSTQDICRICHCEG
 DEESPLITPCRCTGTLRFVHQSLHQWIKSSDTRCCCELCKYDFIMETKPKLRKWEKLMQTTSERRKIFC
 SVTFHVIAITCVVWSLYVLIDRTAEEIKQGNNDNGVLEWPFWTKLVVVAIGFTGGLVFMVYQCKVYVQLWR
 RLKAYNRVIFVQNCPTAKKLEKNFSCNVNTDIKDAVVVVPVPTGANS LPSAEGGPPEVSV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8048_e03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_017923

ORF Size: 816 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_017923.4](#)

RefSeq Size: 1790 bp

RefSeq ORF: 819 bp

Locus ID: 55016

UniProt ID: [Q8TCQ1](#)

Cytogenetics: 4q32.2-q32.3


Domains: RING

Protein Families: Druggable Genome, Transmembrane

MW: 30.7 kDa

Gene Summary: MARCH1 is a member of the MARCH family of membrane-bound E3 ubiquitin ligases (EC 6.3.2.19). MARCH proteins add ubiquitin (see MIM 191339) to target lysines in substrate proteins, thereby signaling their vesicular transport between membrane compartments. MARCH1 downregulates the surface expression of major histocompatibility complex (MHC) class II molecules (see MIM 142880) and other glycoproteins by directing them to the late endosomal/lysosomal compartment (Bartee et al., 2004 [PubMed 14722266]; Thibodeau et al., 2008 [PubMed 18389477]; De Gassart et al., 2008 [PubMed 18305173]). [supplied by OMIM, Mar 2010]

Product images:



Circular map for RC222000