

Product datasheet for **RC222024A1V**

Human PLP2 (NM_002668) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PLP2 (NM_002668) AAV Particle
Tag: Myc-DDK
Symbol: PLP2
Synonyms: A4; A4LSB
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC222024 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGATTCTGAGCGCCTCTCGGCTCCTGGCTGCTGGGCCGCTGCACCAACTTCTCGCGCACTCGAA
AGGGAATCCTCCTGTTTGCTGAGATTATATTATGCCTGGTGATCCTGATCTGCTTCAGTGCCTCCACACC
AGGCTACTCCTCCCTGTCGGTGATTGAGATGATCCTTGCTGCTATTTCTTTGTTGTCTACATGTGTGAC
CTGCACACCAAGATAACCATTATCAACTGGCCCTGGAGTGATTCTTCCGAACCCCTATAGCGGCAATCC
TCTACCTGATCACCTCCATTGTTGCTCTTGTGAGAGAGGAAACCACTCCAAAATCGTCGCAGGGGTACT
GGGCCTAATCGCTACGTGCCTCTTTGGCTATGATGCCTATGTCACCTTCCCGTTCCGGCAGCCAAGACAT
ACAGCAGCCCCCACTGACCCCGCAGATGGCCCGGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC222024 protein sequence
Red=Cloning site Green=Tags(s)

MADSERLSAPGCWAACTNFSRTRKGILLFAEIIILCLVILICFSASTPGYSSLSVIEMILAAIFFVYVYMD
LHTKIPFINWPWSDFFRTLIAAILYLITSIVLVERGNHSHKIVAGVLGLIATCLFGYDAYVTFPVRQPRH
TAAPTDPADGPV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_002668
ORF Size:	456 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_002668.1</u>
RefSeq Size:	1049 bp
RefSeq ORF:	459 bp
Locus ID:	5355
UniProt ID:	<u>Q04941</u>
Cytogenetics:	Xp11.23
MW:	16.7 kDa