

Product datasheet for **RC222034L4V**

RAVER1 (NM_133452) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RAVER1 (NM_133452) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RAVER1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133452
ORF Size:	2268 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222034).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133452.2 , NP_597709.2
RefSeq Size:	3606 bp
RefSeq ORF:	2220 bp
Locus ID:	125950
UniProt ID:	Q8IY67
Cytogenetics:	19p13.2
MW:	79.4 kDa
Gene Summary:	Cooperates with PTBP1 to modulate regulated alternative splicing events. Promotes exon skipping. Cooperates with PTBP1 to modulate switching between mutually exclusive exons during maturation of the TPM1 pre-mRNA (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)