

## Product datasheet for **RC222058L4V**

### SLC34A3 (NM\_080877) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC34A3 (NM_080877) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC34A3
Synonyms:	HHRH; NPTIIc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080877
ORF Size:	1797 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222058).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_080877.1</a> , <a href="#">NP_543153.1</a>
RefSeq Size:	2124 bp
RefSeq ORF:	1800 bp
Locus ID:	142680
UniProt ID:	<a href="#">Q8N130</a>
Cytogenetics:	9q34.3
Protein Families:	Transmembrane
MW:	63.3 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a member of SLC34A transporter family of proteins, and is expressed primarily in the kidney. It is involved in transporting phosphate into cells via sodium cotransport in the renal brush border membrane, and contributes to the maintenance of inorganic phosphate concentration in the kidney. Mutations in this gene are associated with hereditary hypophosphatemic rickets with hypercalciuria. Alternatively spliced transcript variants varying in the 5' UTR have been found for this gene.[provided by RefSeq, Apr 2010]