

Product datasheet for **RC222086L3V**

FOXD2 (NM_004474) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FOXD2 (NM_004474) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FOXD2
Synonyms:	FKHL17; FREAC-9; FREAC9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004474
ORF Size:	1485 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222086).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004474.2 , NP_004465.2
RefSeq Size:	4689 bp
RefSeq ORF:	1488 bp
Locus ID:	2306
UniProt ID:	O60548
Cytogenetics:	1p33
Protein Families:	Druggable Genome, Transcription Factors
MW:	48.6 kDa



[View online »](#)

Gene Summary:

This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined.
[provided by RefSeq, Jul 2008]