

## Product datasheet for **RC222093L4V**

### MTUS1 (NM\_001001925) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	MTUS1 (NM_001001925) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MTUS1
Synonyms:	ATBP; ATIP; ATIP3; ICIS; MP44; MTSG1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001001925
ORF Size:	3648 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222093).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001001925.2</a> , <a href="#">NP_001001925.1</a>
RefSeq Size:	6273 bp
RefSeq ORF:	3651 bp
Locus ID:	57509
UniProt ID:	<a href="#">Q9ULD2</a>
Cytogenetics:	8p22
Protein Families:	Druggable Genome, Transcription Factors, Transmembrane
MW:	135.7 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein which contains a C-terminal domain able to interact with the angiotension II (AT2) receptor and a large coiled-coil region allowing dimerization. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. One of the transcript variants has been shown to encode a mitochondrial protein that acts as a tumor suppressor and participates in AT2 signaling pathways. Other variants may encode nuclear or transmembrane proteins but it has not been determined whether they also participate in AT2 signaling pathways. [provided by RefSeq, Jul 2008]