

## Product datasheet for **RC222145L3V**

### **KCNH6 (NM\_173092) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	KCNH6 (NM_173092) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KCNH6
Synonyms:	ERG-2; ERG2; hERG-2; HERG2; Kv11.2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_173092
ORF Size:	2715 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222145).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_173092.1</a>
RefSeq Size:	3017 bp
RefSeq ORF:	2718 bp
Locus ID:	81033
UniProt ID:	<a href="#">Q9H252</a>
Cytogenetics:	17q23.3
Protein Families:	Druggable Genome, Ion Channels: Other, Transmembrane
MW:	100.3 kDa



[View online »](#)

**Gene Summary:**

Voltage-gated potassium (Kv) channels represent the most complex class of voltage-gated ion channels from both functional and structural standpoints. Their diverse functions include regulating neurotransmitter release, heart rate, insulin secretion, neuronal excitability, epithelial electrolyte transport, smooth muscle contraction, and cell volume. This gene encodes a member of the potassium channel, voltage-gated, subfamily H. This member is a pore-forming (alpha) subunit. Alternative splicing results in multiple transcript variants that encode different isoforms. [provided by RefSeq, Jul 2013]