

Product datasheet for **RC222174L1V**

EXD3 (NM_017820) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EXD3 (NM_017820) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EXD3
Synonyms:	mut-7; Nbr
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_017820
ORF Size:	2628 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222174).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017820.3
RefSeq Size:	2887 bp
RefSeq ORF:	2631 bp
Locus ID:	54932
UniProt ID:	Q8N9H8
Cytogenetics:	9q34.3
MW:	96.4 kDa


[View online »](#)

Gene Summary: Possesses 3'-5' exoribonuclease activity. Required for 3'-end trimming of AGO1-bound miRNAs (By similarity).[UniProtKB/Swiss-Prot Function]