

Product datasheet for **RC222182L2V**

ZNF161 (VEZF1) (NM_007146) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF161 (VEZF1) (NM_007146) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF161
Synonyms:	DB1; ZNF161
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_007146
ORF Size:	1563 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222182).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007146.2
RefSeq Size:	4652 bp
RefSeq ORF:	1566 bp
Locus ID:	7716
UniProt ID:	Q14119
Cytogenetics:	17q22
Domains:	zf-C2H2
Protein Families:	Transcription Factors


[View online »](#)

MW: 56.8 kDa

Gene Summary: Transcriptional regulatory proteins containing tandemly repeated zinc finger domains are thought to be involved in both normal and abnormal cellular proliferation and differentiation. ZNF161 is a C2H2-type zinc finger protein (Koyano-Nakagawa et al., 1994 [PubMed 8035792]). See MIM 603971 for general information on zinc finger proteins.[supplied by OMIM, Sep 2008]