

Product datasheet for **RC222183**

SETDB2 (NM_031915) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SETDB2 (NM_031915) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SETDB2
Synonyms:	C13orf4; CLLD8; CLLL8; KMT1F
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC222183 representing NM_031915
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGGGAGAAAAAATGGCGATGCAAAAACCTTTCTGGATGGAGCTAGAAGATGATGGAAAAGTGGACTTCA
 TTTTGAACAAGTACAAAATGTGCTGCAGTCACTGAAACAAAAGATCAAAGATGGGTCTGCCACCAATAA
 AGAATACATCCAAGCAATGATTCTAGTGAATGAAGCAACTATAATTAACAGTTCAACATCAATAAAGGGA
 GCATCACAGAAAGAAGTGAATGCCCAAAGCAGTGATCCTATGCCTGTGACTCAGAAGGAACAGGAAAAACA
 AATCCAATGCATTTCCCTCTACATCATGTGAAAACCTCTTCCAGAAGACTGTACATTTCTAACAACAGG
 AAATAAGGAAATCTCTCTCTTGAAGATAAAGTTGTAGACTTTAGAGAAAAAGACTCATCTTCGAATTTA
 TCTTACCAAAGTCATGACTGCTCTGGTGTCTGTGATGAAAATGCCACTGAACTGAAGGGAGAAAAACC
 CTCTGCAGCTGCCAATCAAATGTCACTTCAAAGACGACATGCAAAGACAACTCTCATTCTTCAGCACT
 CCACGTGAGTTATAAAACCCCTTGTGGAAGGAGTCTACGAAACGTGGAGGAAGTTTTTCGTTACCTGCTT
 GAGACAGAGTGTAACTTTTATTACAGATAACTTTTTCTTCAATACCTATGTTCAAGTTGGCTCGGAATT
 ACCCAAAGCAAAAAGAAGTTGTTTCTGATGTGGATATTAGCAATGGAGTGAATCAGTGCCCATTTCTTT
 CTGTAATGAAATTGACAGTAGAAAGCTCCACAGTTTAAGTACAGAAAGACTGTGTGGCCTCGAGCATAT
 AATCTAACCAACTTTTCCAGCATGTTTACTGATTCCTGTGACTGCTCTGAGGGCTGCATAGACATAACAA
 AATGTGCATGTCTTCACTGACAGCAAGGAATGCCAAAACCTCCCTTGTCAAGTGACAAAATAACAC
 TGGATATAAATATAAAGACTACAGAGACAGATTCTACTGGCATTATGAATGCAGCCTTTTGTGCAAA
 TGTAAATCGACAATTGTGTCAAAACCGAGTTGTCCAACATGGTCTCAAGTGAGGTTACAGGTGTTCAAAA
 CTGAGCAGAAGGGATGGGTGTACGCTGTCTAGATGACATTGACAGAGGGACATTTGTTTGCAATTTATTC
 AGGAAGATTACTAAGCAGAGCTAACACTGAAAAATCTTATGGTATTGATGAAAACGGGAGAGATGAGAAT
 ACTATGAAAAATATATTTTCAAAAAAGAGGAAATTAGAAGTTGCATGTTCAAGATTGTGAAGTTGAAGTTC
 TCCCATTAGGATTGGAACACATCCTAGAACTGCTAAACTGAGAAATGTCCACCAAAGTTCAGTAATAA
 TCCCAAGGAGCTTACTGTGGAACGAAATATGATAATTTCAAGAATTCAATATCATTACAGTTATTAGA
 GATCCTGAATCCAAGACAGCCATTTTCAACACAATGGGAAAAAATGGAATTTGTTTCTCGGAGTCTG
 TCACTCCAGAAGATAATGATGGATTTAAACCACCCGAGAGCATCTGAACTCTAAACCAAGGGAGCACA
 AAAGGACTCAAGTTCAAACCATGTTGATGAGTTTGAAGATAATCTGCTGATTGAATCAGATGTGATAGAT
 ATAATAAATATAGAGAAGAACTCCACCAAGGAGCAGATGTAACCAGGCGACCACATTGGATAATCAGA
 ATATTAATAAAGGCAATTGAGGTTCAAATTCAGAAACCCCAAGAGGGACGATCTACAGCATGTCAAAGACA
 GCAGGTATTTTGTGATGAAGAGTTGCTAAGTGAACCAAGAATACTTCATCTGATTCTCTAACAAGTTC
 AATAAAGGGAATGTGTTTTATTGGATGCCACAAAAGAAGGAAATGTGCGCCGCTTCCTTAATCATAGTT
 GTTGCCCAAACTCTTGGTACAGAATGTTTTGTAGAAACACACAACAGGAATTTCCATTGGTGGCATT
 CTTACCAACAGGTATGTGAAAGCAAGAACAGAGCTAACATGGGATTATGGCTATGAAGCTGGGACTGTG
 CCTGAGAAGGAAATCTTCTGCCAATGTGGGGTTAATAAATGTAGAAAAAATATTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC222183 representing NM_031915
 Red=Cloning site Green=Tags(s)


MGEKNGDAKTFWMELEDDGKVDIFEFQVQNVLSLQSLKQIKDGSATNKEYIQAMILVNEATIINSSTSIK
 ASQKEVNAQSSDPMPTQKEQENKSNAPSTSCENSFPEDCTFLTGTNKEILSLEDKVVDFREKSSNL
 SYQSHDCSGACLMKMLNLKGENPLQLPIKCHFQRRHAKTNSHSSALHVSYKTPCGRSLRNVEEVFRYLL
 ETECNFLFTDNFSFNTYVQLARNYPKQKEVSDVDISNGVESVPISFCNEIDSRKLPQFKYRKTVPWAY
 NLTNFSSMFTDSCDCSEGCIDITKCAQLTARNAKTSPLSSDKITTGKYKRLQRQIPTGIYECSLCK
 CNRQLCQNRVVQHGPQVRLQVFKTEQKGWVRCCLDDIDRGTFCIYSGRLLSRANTEKSYGIDENGRDEN
 TMKNIFSKKRKLEVACSDCEVEVLPLGLETHPRTAKTEKCPPKFSNNPKELTVETKYDNISRIQYHSVIR
 DPESKTAIFQHNGKKMEFVSSESVTPEDNDGFKPPREHLNSKTGAQKSSSNHVDEFEDNLLIESDVID
 ITKYREETPPSRCNQATTLNQNIIKKAIEVQIQKPQEGRSTACQRRQVFCDEELLSETKNTSSDSLTKF
 NKGNVFLLDATKEGNVGRFLNHSCCPNLLVQNVFVETHNRNPLVAFFTNRYVKARTELTDYGYEAGTV
 PEKEIFCQCGVNKCRKKIL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mg3548_e01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_031915


ORF Size: 2157 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_031915.3</u>
RefSeq Size:	3307 bp
RefSeq ORF:	2160 bp
Locus ID:	83852
UniProt ID:	<u>Q96T68</u>
Cytogenetics:	13q14.2
Domains:	SET, MBD, PreSET, Pre-SET
Protein Families:	Druggable Genome
Protein Pathways:	Lysine degradation
MW:	81.7 kDa
Gene Summary:	This gene encodes a member of a family of proteins that contain a methyl-CpG-binding domain (MBD) and a SET domain and function as histone methyltransferases. This protein is recruited to heterochromatin and plays a role in the regulation of chromosome segregation. This region is commonly deleted in chronic lymphocytic leukemia. Naturally-occurring readthrough transcription occurs from this gene to the downstream PHF11 (PHD finger protein 11) gene. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016]

Product images:



Circular map for RC222183