

Product datasheet for **RC222213A1V**

Human OR51B2 (NM_033180) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human OR51B2 (NM_033180) AAV Particle
Tag:	Myc-DDK
Symbol:	OR51B2
Synonyms:	HOR5'Beta3; OR51B1P
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC222213 representing NM_033180 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTGGCCCAATATTACTGCAGCCCCTTTTGTGCTGACTGGTTTTCCAGGGCTGGAGGCAGCTCATCACT
GGATCTCCATCCCCTTCTTTGCTGTTTATGTGTGCATCCTTCTGGGCAATGGCATGCTCCTCTACCTCAT
CAAGCATGACCACAGTCTTCATGAGCCCATGTACTACTTCCTCACCATGCTGGCAGGCACAGACCTCATG
GTGACATTGACCACGATGCCTACTGTAATGGGCATCCTATGGTGAATCACAGGGAGATTAGCAGTGTGG
GCTGCTTCCTACAGGCTTACTTTATTCACCTCCCTTTCTGTTGTGGAATCAGGTTCCCTCCTGGCAATGGC
ATATGATCGTTTCATTGCCATCCGCAATCCTTTGAGATATGCTTCCATTCTACCAATACTAGAGTCATA
GCGTTAGGAGTGGGAGTGTCTAAGGGGTTTTGTATCCATCCTGCCTGTAATTTTTCGCTCTTTTTTCAT
TTTCATATTGCAAATCTCATGTTATCACACGTGCTTTCTGCCTCCACCAAGAAATCATGAGACTGGCTTG
TGCTGACATAACTTTCAATAGACTTTACCCTGTAATTTTGTATCTCTTTAACAATCTTCTAGACTCTCTG
ATCATCCTCTTCTCTATATTCTAATTCTTAATACTGTCATAGGCATTGCCTCTGGTGAAGAGAGAGCCA
AAGCCCTCAATACCTGTATCTCCACATTAGTTGTGTTCTTATCTTCTATGTTACGGTGTATGGGTTTGAC
ATTCATTTACAGATTTGGGAAGAATGTGCCAGAGGTTGTCCACATTATCATGAGTTACATCTACTTCCTC
TTTCCTTCTTTAATGAACCCTGTCATCTACAGCATCAAAACCAAGCAAAATACAATATGGCATTATCCGCC
TTTTATCTAAACATAGGTTTAGTAGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC222213 representing NM_033180 Red=Cloning site Green=Tags(s) MWPNITAAPFLLTGFPGLEAAHHWISIPFFAVYVCILLGNGMLLYLIKHDHSLHEPMYYFLTMLAGTDLM VLTMTPTVMGILWVNHREISSVGCFLQAYFIHSLSVVESGSLAMAYDRFIAIRNPLRYASILTNTRVI ALGVGVFLRGFVSILPVILRLFSFSYCKSHVITRAFCLHQEIMRLACADITFNRLYPVILISLTIFLDSL IILFSYILILNTVIGIASGEERAKALNTCISHISCVLIFYVTVMGLTFIYRFGKNVPEVVHIIMSYIYFL FPSLMNPVIYSIKTKQIQYGIIRLLSKHRFSS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_033180
ORF Size:	936 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_033180.1</u>
RefSeq Size:	939 bp
RefSeq ORF:	939 bp
Locus ID:	79345
UniProt ID:	<u>Q9Y5P1</u>
Cytogenetics:	11p15.4
MW:	35.2 kDa