

Product datasheet for **RC222220L3V**

DNAJC14 (NM_032364) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DNAJC14 (NM_032364) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DNAJC14
Synonyms:	DNAJ; DRIP78; HDJ3; LIP6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032364
ORF Size:	2106 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222220).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032364.5 , NP_115740.5
RefSeq Size:	3330 bp
RefSeq ORF:	2109 bp
Locus ID:	85406
UniProt ID:	Q6Y2X3
Cytogenetics:	12q13.2
Domains:	Dnaj
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

MW: 78.6 kDa

Gene Summary: Regulates the export of target proteins, such as DRD1, from the endoplasmic reticulum to the cell surface.[UniProtKB/Swiss-Prot Function]