

## Product datasheet for **RC222233L3V**

### AGO3 (NM\_024852) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | AGO3 (NM_024852) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | AGO3   |
| Synonyms:                 | EIF2C3   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_024852  |
| ORF Size:                 | 2580 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC222233).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_024852.2</a>  |
| RefSeq Size:              | 3578 bp  |
| RefSeq ORF:               | 2583 bp  |
| Locus ID:                 | 192669   |
| UniProt ID:               | <a href="#">Q9H9G7</a>   |
| Cytogenetics:             | 1p34.3   |
| Domains:                  | PAZ, Piwi  |
| Protein Families:         | Druggable Genome   |



[View online »](#)

**MW:** 97.2 kDa

**Gene Summary:** This gene encodes a member of the Argonaute family of proteins which play a role in RNA interference. The encoded protein is highly basic, contains a PAZ domain and a PIWI domain, and may play a role in short-interfering-RNA-mediated gene silencing. This gene is located on chromosome 1 in a tandem cluster of closely related family members including argonaute 4 and eukaryotic translation initiation factor 2C, 1. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]