

Product datasheet for **RC222240A1V**

Human PACRG (NM_001080378) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PACRG (NM_001080378) AAV Particle
Tag:	Myc-DDK
Symbol:	PACRG
Synonyms:	GLUP; HAK005771; PACRG2.1; PARK2CRG
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC222240 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTGGCAGAAAAAGAGACCCTGAGCTTAAACAAATGCCAGACAAGATGCCGAAGAGGACCAAGCTGC
TGGCACAACAGCCGCTCCCGGTGCACGAGCCTCACTCTCTGGTTTCTGAGGGTTTCACAGTCAAAGCCAT
GATGAAAACTCAGTCGTGAGAGGCCCTCCAGCTGCAGGGCATTAAAGAAAGACCAACCAAGCCACACA
GCATTTGAAAAATTCTATGAGCGAGGTGACTTCCCAATTGCCCTTGAGCATGATTCGAAAGGAAACAAAA
TCGCCTGGAAGGTAGAAATTGAGAAGCTGGATTACCATCATTATCTGCCTCTGTTTTTGGTGGCTTTG
TGAAATGACATTTCCCTATGAGTTTTTTGCTCGGCAAGGAATCCACGACATGCTGGAACACGGTGGGAAC
AAGATCCTACCTGTCCTTCCACAGCTCATTATCCCGATAAAAAATGCCTTGAACCTCCGAAACCGACAGG
TCATCTGTGCTCACTCTCAAGGTCCTCCAGCATCTGGTTGTGTCAGCTGAGATGGTGGCAAGGCCTTGGT
GCCTTATTACCGTCAAATCCTCCTGCTGAACATCTTTAAGAATATGAATGTGAACCTCCGAGACGGC
ATTGACTACAGCCAGCAGAAGAGGGAGAACATTGGGGACTTGATCCAGGAGACACTGGAGGCCTTCGAGC
GCTACGGAGGAGAAAAATGCCTTTATCAACATTAAGTACGTGGTCCCAACCTACGAGTCTTGCTTGCTAAA
C

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC222240 protein sequence Red=Cloning site Green=Tags(s) MVAEKETLSLNKCPDKMPKRTKLLAQQLPVHQPHSLVSEGFTVKAMMKNSVVRGPPAAGAFKERPTKPT AFRKFYERGDFPIALEHDSKGNKIAWKVEIEKLDYHHYLPLFFDGLCEMTFPYEFFARQGIHDMLEHGGN KILPVLPLQLIPIKNALNLRNRQVICVTLKVLQHLVSAEMVGKALVPYYRQILPVLNIFKNMNVNSGDG IDYSQQKRENIGDLIQETLEAFERYGGENAFINIKYVVPTYESCLLN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001080378
ORF Size:	771 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001080378.1, NP_001073847.1</u>
RefSeq Size:	1585 bp
RefSeq ORF:	774 bp
Locus ID:	135138
UniProt ID:	<u>Q96M98</u>
Cytogenetics:	6q26
MW:	29.3 kDa