

Product datasheet for **RC222271L4V**

Hemopexin (HPX) (NM_000613) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hemopexin (HPX) (NM_000613) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Hemopexin
Synonyms:	HX
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000613
ORF Size:	1386 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222271).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000613.1
RefSeq Size:	1389 bp
RefSeq ORF:	1389 bp
Locus ID:	3263
UniProt ID:	P02790
Cytogenetics:	11p15.4
Domains:	hemopexin
Protein Families:	Secreted Protein



[View online »](#)

MW: 51.68 kDa

Gene Summary: This gene encodes a plasma glycoprotein that binds heme with high affinity. The encoded protein is an acute phase protein that transports heme from the plasma to the liver and may be involved in protecting cells from oxidative stress. [provided by RefSeq, Apr 2009]