

Product datasheet for **RC222274**

MPIG6B (NM_025260) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: MPIG6B (NM_025260) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: MPIG6B
Synonyms: C6orf25; G6b; G6b-B; NG31; THAMY
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC222274 representing NM_025260
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTGTGTTTCTGCAGCTGCTACCGCTGCTGCTCTCGAGGGCCCAAGGGAACCCTGGGGCTTCTCTGG
 ACGGCCGCCCTGGGGACCGGGTGAATCTCTCTCGCGAGGAGTCTCTCATCCCATCCGCTGGGTCTGGGC
 ACCCAGCTTCCCGCCTGCAAGGGCCTGTCAAAGGACGCCGACCGATCCTGTGGCCTCTTCGAGCGGG
 ACCCCACCGTGCCTCCCTCCAGCCTTTCGTGCGCCGCTACGCTCCCTGGACTCTGGTATCCGGCGGC
 TGGAGCTCCTCTTGAAGCGGGGGACTCGGGCACTTTTTTCTGCAAGGGCCGCCACGAGGACGAGAGCCG
 TACAGTGTTCACGTGCTGGGGACAGGACCTATTGCAAGGCCCCCGGGCCTACCCATGGGTCCGTGTAT
 CCCAGCTCCTGATCCCGCTGCTGGGCGCTGGGTTGGTGCTCGGACTGGGAGCTTTGGGCCTGGTCTGGT
 GGCTGCACAGGCGCTGCCCCGCAACCGATTGACCACTCCCTAGATTTGCTCTGTCCCCCCACATAG
 CTCCACTTGTGAAAACCGAGCCCCAGAGGCCAGTAAAGGAGGAAGAGCCCAAGATTCCAGGGGACCTGGA
 CCAGGAACCGAGCCTGCTCTATGCGGATCTGGACCATCTAGCCCTCAGCAGGCCCCCGCGGCTGTCCACA
 GCGGACCCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RC222274 representing NM_025260
Red=Cloning site Green=Tags(s)

MAVFLQLLLPLLLSRAQGNPGASLDGRPGDRVNLSCGGVSHPIRWVWAPSFPAKGLSKGRRPILWASSSG
 TPTVPPLQPFVGRRLSLDSGIRRELELLSAGDSGTFCKGRHEDESRTLHLVLDRTYCKAPGPTHGSVY
 PQLLIPLLGAGLVGLGALGLVWWLHRRLLPPQPIRPLPRFALSPHSSTCENRAPEASKGGRAQDSRGP
 PGTEPALCGSGPSSPQQAPPAVHSGPC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_025260

ORF Size: 711 bp


OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_025260.4
RefSeq Size:	2435 bp
RefSeq ORF:	714 bp
Locus ID:	80739
UniProt ID:	O95866
Cytogenetics:	6p21.33
Protein Families:	Druggable Genome, Transmembrane
MW:	25 kDa
Gene Summary:	This gene is a member of the immunoglobulin (Ig) superfamily and is located in the major histocompatibility complex (MHC) class III region. The protein encoded by this gene is a glycosylated, plasma membrane-bound cell surface receptor, but soluble isoforms encoded by some transcript variants have been found in the endoplasmic reticulum and Golgi before being secreted. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC222274