

## Product datasheet for **RC222280L4V**

### ADAM12 (NM\_021641) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ADAM12 (NM_021641) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADAM12
Synonyms:	ADAM12-OT1; CAR10; MCMP; MCMPMItna; MLTN; MLTNA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021641
ORF Size:	2214 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222280).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021641.3</a>
RefSeq Size:	3355 bp
RefSeq ORF:	2217 bp
Locus ID:	8038
UniProt ID:	<a href="#">O43184</a>
Cytogenetics:	10q26.2
Protein Families:	Druggable Genome, Protease, Secreted Protein, Transmembrane
MW:	77.6 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of a family of proteins that are structurally related to snake venom disintegrins and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. Expression of this gene has been used as a maternal serum marker for pre-natal development. Alternative splicing results in multiple transcript variants encoding different isoforms. Shorter isoforms are secreted, while longer isoforms are membrane-bound form. [provided by RefSeq, Jan 2014]