

Product datasheet for **RC222281L2V**

Raptor (RPTOR) (NM_020761) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Raptor (RPTOR) (NM_020761) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Raptor
Synonyms:	KOG1; Mip1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_020761
ORF Size:	4005 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222281).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020761.1
RefSeq Size:	6431 bp
RefSeq ORF:	4008 bp
Locus ID:	57521
UniProt ID:	Q8N122
Cytogenetics:	17q25.3
Domains:	WD40
Protein Pathways:	Insulin signaling pathway, mTOR signaling pathway



[View online »](#)

MW: 148.9 kDa

Gene Summary: This gene encodes a component of a signaling pathway that regulates cell growth in response to nutrient and insulin levels. The encoded protein forms a stoichiometric complex with the mTOR kinase, and also associates with eukaryotic initiation factor 4E-binding protein-1 and ribosomal protein S6 kinase. The protein positively regulates the downstream effector ribosomal protein S6 kinase, and negatively regulates the mTOR kinase. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]