

Product datasheet for **RC222282L4V**

KRT26 (NM_181539) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KRT26 (NM_181539) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KRT26
Synonyms:	CK26; K25B; K25IRS2; K26; KRT25B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_181539
ORF Size:	1404 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222282).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_181539.3 , NP_853517.2
RefSeq Size:	1748 bp
RefSeq ORF:	1407 bp
Locus ID:	353288
UniProt ID:	Q7Z3Y9
Cytogenetics:	17q21.2
MW:	51.9 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is a member of the keratin superfamily. This keratin protein is a type I keratin that is specific for the inner root sheath of the hair follicle. This gene exists in a cluster with other keratin genes on chromosome 17q21. [provided by RefSeq, Jul 2009]