

## Product datasheet for **RC222287L3V**

### MTMR3 (NM\_021090) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	MTMR3 (NM_021090) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MTMR3
Synonyms:	FYVE-DSP1; ZFYVE10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_021090
ORF Size:	3594 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222287).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021090.2</a>
RefSeq Size:	5963 bp
RefSeq ORF:	3597 bp
Locus ID:	8897
UniProt ID:	<a href="#">Q13615</a>
Cytogenetics:	22q12.2
Domains:	FYVE, PTPc_motif
Protein Families:	Druggable Genome, Phosphatase



[View online »](#)

**MW:** 133.4 kDa

**Gene Summary:** This gene encodes a member of the myotubularin dual specificity protein phosphatase gene family. The encoded protein is structurally similar to myotubularin but in addition contains a FYVE domain and an N-terminal PH-GRAM domain. The protein can self-associate and also form heteromers with another myotubularin related protein. The protein binds to phosphoinositide lipids through the PH-GRAM domain, and can hydrolyze phosphatidylinositol(3)-phosphate and phosphatidylinositol(3,5)-biphosphate in vitro. The encoded protein has been observed to have a perinuclear, possibly membrane-bound, distribution in cells, but it has also been found free in the cytoplasm. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]