

## Product datasheet for **RC222307L2V**

### **FKBP9 (NM\_007270) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	FKBP9 (NM_007270) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FKBP9
Synonyms:	FKBP60; FKBP63; PPIase
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_007270
ORF Size:	1710 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222307).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007270.2</a>
RefSeq Size:	3425 bp
RefSeq ORF:	1713 bp
Locus ID:	11328
UniProt ID:	<a href="#">O95302</a>
Cytogenetics:	7p14.3
Domains:	EFh
Protein Families:	Druggable Genome


[View online »](#)

**MW:** 62.9 kDa

**Gene Summary:** PPlases accelerate the folding of proteins during protein synthesis.[UniProtKB/Swiss-Prot Function]