

Product datasheet for **RC222318L4V**

RPS4Y2 (NM_001039567) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RPS4Y2 (NM_001039567) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RPS4Y2
Synonyms:	RPS4Y2P
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001039567
ORF Size:	789 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222318).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001039567.1
RefSeq Size:	792 bp
RefSeq ORF:	792 bp
Locus ID:	140032
UniProt ID:	Q8TD47
Cytogenetics:	Yq11.223
MW:	29.1 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is a ribosomal protein that is highly similar to RPS4Y1. This gene is located in the male-specific region of the Y chromosome. [provided by RefSeq, Aug 2012]