

Product datasheet for **RC222328L3V**

NAALADL2 (NM_207015) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NAALADL2 (NM_207015) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NAALADL2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_207015
ORF Size:	2385 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222328).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_207015.1 , NP_996898.1
RefSeq Size:	4928 bp
RefSeq ORF:	2388 bp
Locus ID:	254827
UniProt ID:	Q58DX5
Cytogenetics:	3q26.31
Protein Families:	Protease, Transmembrane
MW:	88.5 kDa
Gene Summary:	May be catalytically inactive.[UniProtKB/Swiss-Prot Function]



[View online »](#)