

## Product datasheet for **RC222329L4V**

### EXOC7 (NM\_015219) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	EXOC7 (NM_015219) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EXOC7
Synonyms:	2-5-3p; BLOM4; EX070; EXO70; Exo70p; EXOC1; NEDSEBA; YJL085W
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015219
ORF Size:	1959 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222329).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015219.2</a>
RefSeq Size:	4687 bp
RefSeq ORF:	1962 bp
Locus ID:	23265
UniProt ID:	<a href="#">Q9UPT5</a>
Cytogenetics:	17q25.1
Domains:	Exo70
Protein Families:	Druggable Genome


[View online »](#)

**Protein Pathways:** Insulin signaling pathway

**MW:** 74.5 kDa

**Gene Summary:** The protein encoded by this gene is a component of the exocyst complex. The exocyst complex plays a critical role in vesicular trafficking and the secretory pathway by targeting post-Golgi vesicles to the plasma membrane. The encoded protein is required for assembly of the exocyst complex and docking of the complex to the plasma membrane. The encoded protein may also play a role in pre-mRNA splicing through interactions with pre-mRNA-processing factor 19. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 4. [provided by RefSeq, Nov 2011]