

Human TAS2R3 (NM 016943) AAV Particle

Product Type:	AAV Particles
Product Name:	Human TAS2R3 (NM_016943) AAV Particle
Tag:	Myc-DDK
Symbol:	TAS2R3
Synonyms:	T2R3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC222347 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

[illegible]

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

[View online »](#)

Protein Sequence:	>RC222347 protein sequence Red=Cloning site Green=Tags(s) MMGLTEGVFLILSGTQFTLGILVNCFIELVNGSSWFKTKRMSLSDFIITTLALLRIILLCIILTDSFLIE FSPNTHDSGIIMQIIDVSWFTFNHLSIWLATCLGVLYCLKIASFSHPTFLWLKWRVSRVMVMMLLGALLL SCGSTASLINEFKLYSVFRGIEATRNVEHFRKKRSEYYLIHVLGTLWYLPPLIVSLASYSLLIFSLGRH TRQMLQNGTSSRDPTEAHKRAIRIILSFFFLFLLYFLAFLIASFGNFLPKTKMAKMIGEVMTFYPAGH SFILILGNSKLKQTFVVMRLRCESGHLKPGSGKPIFS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_016943
ORF Size:	948 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_016943.2 , NP_058639.1
RefSeq Size:	1101 bp
RefSeq ORF:	951 bp
Locus ID:	50831
UniProt ID:	Q9NYW6
Cytogenetics:	7q34
MW:	35.9 kDa