

Product datasheet for **RC222349L4V**

TNS4 (NM_032865) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TNS4 (NM_032865) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TNS4
Synonyms:	CTEN; PP14434
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032865
ORF Size:	2145 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222349).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032865.5
RefSeq Size:	4090 bp
RefSeq ORF:	2148 bp
Locus ID:	84951
UniProt ID:	Q8IZW8
Cytogenetics:	17q21.2
Domains:	SH2, PID
MW:	76.6 kDa



[View online »](#)

Gene Summary:

May be involved in cell migration, cartilage development and in linking signal transduction pathways to the cytoskeleton (By similarity). May promote apoptosis, via its cleavage by caspase-3.[UniProtKB/Swiss-Prot Function]