

### **Product datasheet for RC222399L1**

# DLL1 (NM\_005618) Human Tagged Lenti ORF Clone

#### **Product data:**

**Product Type:** Expression Plasmids

**Product Name:** DLL1 (NM\_005618) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: DLL1

Synonyms: Delta; DELTA1; DL1; NEDBAS

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC222399).

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.

**ACCN:** NM\_005618

ORF Size: 2169 bp



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



#### DLL1 (NM\_005618) Human Tagged Lenti ORF Clone - RC222399L1

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 005618.3</u>

 RefSeq Size:
 3366 bp

 RefSeq ORF:
 2172 bp

 Locus ID:
 28514

 UniProt ID:
 000548

 Cytogenetics:
 6q27

**Domains:** DSL, EGF CA, EGF, EGF

Protein Families: Adult stem cells, Cancer stem cells, ES Cell Differentiation/IPS, Stem cell relevant signaling -

DSL/Notch pathway, Transmembrane

**Protein Pathways:** Notch signaling pathway

**MW:** 77.9 kDa


**Gene Summary:** DLL1 is a human homolog of the Notch Delta ligand and is a member of the

delta/serrate/jagged family. It plays a role in mediating cell fate decisions during


hematopoiesis. It may play a role in cell-to-cell communication. [provided by RefSeq, Jul 2008]



## **Product images:**



Circular map for RC222399L1



Double digestion of RC222399L1 using Sgfl and Mlul