

## Product datasheet for **RC222411L4V**

### **MGEA5 (OGA) (NM\_012215) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	MGEA5 (OGA) (NM_012215) Human Tagged ORF Clone Lentiviral Particle
Symbol:	OGA
Synonyms:	MEA5; MGEA5; NCOAT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012215
ORF Size:	2748 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222411).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012215.1</a>
RefSeq Size:	5220 bp
RefSeq ORF:	2751 bp
Locus ID:	10724
UniProt ID:	<a href="#">O60502</a>
Cytogenetics:	10q24.32
MW:	102.9 kDa



[View online »](#)

**Gene Summary:**

The dynamic modification of cytoplasmic and nuclear proteins by O-linked N-acetylglucosamine (O-GlcNAc) addition and removal on serine and threonine residues is catalyzed by OGT (MIM 300255), which adds O-GlcNAc, and MGEA5, a glycosidase that removes O-GlcNAc modifications (Gao et al., 2001 [PubMed 11148210]).[supplied by OMIM, Mar 2008]