

Product datasheet for RC222426L1

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

RUNX1T1 (NM_175635) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: RUNX1T1 (NM_175635) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: RUNX1T1

Synonyms: AML1-MTG8; AML1T1; CBFA2T1; CDR; ETO; MTG8; ZMYND2

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC222426).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_175635

ORF Size: 1701 bp





RUNX1T1 (NM_175635) Human Tagged Lenti ORF Clone - RC222426L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 175635.1</u>

RefSeq Size:3233 bpRefSeq ORF:1704 bp

Locus ID: 862

UniProt ID: Q06455

Cytogenetics: 8q21.3

Protein Families: Transcription Factors

Protein Pathways: Acute myeloid leukemia, Pathways in cancer

MW: 63 kDa

Gene Summary: This gene encodes a member of the myeloid translocation gene family which interact with

DNA-bound transcription factors and recruit a range of corepressors to facilitate

transcriptional repression. The t(8;21)(q22;q22) translocation is one of the most frequent karyotypic abnormalities in acute myeloid leukemia. The translocation produces a chimeric gene made up of the 5'-region of the runt-related transcription factor 1 gene fused to the 3'-


region of this gene. The chimeric protein is thought to associate with the nuclear

corepressor/histone deacetylase complex to block hematopoietic differentiation. Alternative


splicing results in multiple transcript variants. [provided by RefSeq, Nov 2010]



Product images:



Circular map for RC222426L1



Double digestion of RC222426L1 using Sgfl and Mlul $\,$