

Product datasheet for **RC222431L3V**

Nebulette (NEBL) (NM_006393) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nebulette (NEBL) (NM_006393) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NEBL
Synonyms:	bA165O3.1; C10orf113; LASP2; LNEBL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006393
ORF Size:	3042 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222431).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006393.1 , NP_006384.1
RefSeq Size:	8034 bp
RefSeq ORF:	3045 bp
Locus ID:	10529
UniProt ID:	O76041
Cytogenetics:	10p12.31
MW:	116.3 kDa



[View online »](#)

Gene Summary:

This gene encodes a nebulin like protein that is abundantly expressed in cardiac muscle. The encoded protein binds actin and interacts with thin filaments and Z-line associated proteins in striated muscle. This protein may be involved in cardiac myofibril assembly. A shorter isoform of this protein termed LIM nebulin is expressed in non-muscle cells and may function as a component of focal adhesion complexes. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010]