

## Product datasheet for **RC222445L3V**

### **C1orf173 (ERICH3) (NM\_001002912) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	C1orf173 (ERICH3) (NM_001002912) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C1orf173
Synonyms:	C1orf173
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001002912
ORF Size:	4590 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222445).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001002912.3</a> , <a href="#">NP_001002912.3</a>
RefSeq Size:	7165 bp
RefSeq ORF:	4593 bp
Locus ID:	127254
UniProt ID:	<a href="#">Q5RHP9</a>
Cytogenetics:	1p31.1
MW:	168.3 kDa


[View online »](#)