

## Product datasheet for **RC222450L4V**

### **ABLIM1 (NM\_001003407) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	ABLIM1 (NM_001003407) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ABLIM1
Synonyms:	ABLIM; abLIM-1; LIMAB1; LIMATIN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001003407
ORF Size:	2154 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222450).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001003407.1</a>
RefSeq Size:	7620 bp
RefSeq ORF:	2157 bp
Locus ID:	3983
UniProt ID:	<a href="#">O14639</a>
Cytogenetics:	10q25.3
Protein Pathways:	Axon guidance
MW:	80.9 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a LIM zinc-binding domain-containing protein that binds to actin filaments and mediates interactions between actin and cytoplasmic targets. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jun 2017]