

Product datasheet for RC222465

OR13D1 (NM_001004484) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR13D1 (NM_001004484) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OR13D1
Synonyms:	OR9-15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC222465 representing NM_001004484 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTACAGATTTACAGATTTTGATGTATCAAACATTTCAATTTACCTGAATCATGTCCTTTTCTATACTA
CCCAGCAGGCAGGTGACCTAGAACACATGGAGACAAGAAATTAATCTGCCATGACTGAATCTTTCTGGT
GGGGCTTTCCCAATATCCAGAGCTCCAGCTTTTCTGTTCTGCTCTGCCTCATCATGTACATGATAATC
CTCCTGGGAAATAGCCTCCTCATTATCATCACCATCTGGATTCTCGCCTCCATACTCCCATGTATTCT
TTCTTGAAACCTCTCATTCTTGACATCTGTTACACATCCTCATCCATTCTCCAATGCTTATTATATT
TATGCTGAGAGAAAATCCATCTCCTTCATTGGCTGTGCTCTGCAGATGGTTGTGTCCCTTGGCTTGGC
TCCACTGAGTGTCTCCTCTGGCTGTGATGGCCTATGACCACTATGTGGCCATCTGCAACCCACTGAGGT
ACTCCATCATCATGAACGGAGTGTGTATGTGCAATGGCTGCATGGTCTGGATCATAGGCTGTCTGAC
CTCCCTATTGCAACAGTTCTGACAATGATGTTGCCTTTCTGTGGGAATAATGTCATTGATCATATTACC
TGTGAAATTTGGCCCTTCTAAACTTGTTTGTTCAGATATCACCATCAATGTGCTTATCATGACAGTGA
CAAATATTGTTTCACTGGTGATTCTTCTACTGTTAATTTTCATCTCCTATGTGTTATTCTCTCTCCAT
CCTGAGAATTAATTGTGCTGAGGGAAGAAAGAAAGCCTTCTACCTGTTACGCGCACTCGATTGTGGTC
ATCTTATTCTACGGTTCAGCCCTTTTATGTACATGAAACCAAGTCAAAGAACAATAACATCATGTG
AGATTATTGGCTGTCTTATGGAGTGGTAAGCCCAATGTTAAATCCCATCATCTATAGCCTCAGGAATAA
AGAGGTCAAAGAGGCTGTAAAGAAAGTCTGAGCAGACATCTGCATTTATTGAAAATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001004484.1](#), [NP_001004484.1](#)

RefSeq Size: 1041 bp

RefSeq ORF: 1041 bp

Locus ID: 286365

UniProt ID: [Q8NGV5](#)

Cytogenetics: 9q31.1


Protein Families: Transmembrane

Protein Pathways: Olfactory transduction

MW: 38.9 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC222465