

Product datasheet for **RC222474A1V**

Human FBXO4 (NM_033484) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human FBXO4 (NM_033484) AAV Particle
Tag: Myc-DDK
Symbol: FBXO4
Synonyms: FBX4
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC222474 representing NM_033484
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGGAAGCGAGCCGCGCAGCGGAACAACTCGCCGCCGCCCTTCAGCGACTGGGGCCGCCTGG
AGGCGGCCATCCTCAGCGGCTGGAAGACCTTCTGGCAGTCAGTGAGCAAGGAGAGGGTGGCGCTACGAC
CTCAGGGGAGGAGGTGGATGAGCGGCCAGCACCTGACGCGGCTGCCGATTGATGTACAGCTATATATT
TTGTCCTTTCTTCACCTCATGATCTGTGTGAGTTGGGAAGTACAAATCATTATTGGAATGAACTGTAA
GAGATCCAATTCTGTGGAGATACTTTTGTGGAGGATCTTCTTCTGGTCTTCTGTTGACTGGAAGTC
TCTTCCAGATCTAGAAATCTTAAAAAGCCTATATCTGAGGTCACTGATGGTGCATTTTTTGACTACATG
GCAGTCTATAGAATGTGCTGTCCATACACAAGAAGAGCTTCAAAATCCAGCCGTCCTATGTATGGAGCTG
TCACTTCTTTTTTACACTCCCTGATCATTGAGATGAACACGATTTGCTATGTTTGGACCAGGTTTGA
AGAATTGAATACCTCTTTGGTGTGAGCTTGTGTCTTCAGAGGAACCTTGCCTAACAGCTGGTTTGCCT
CAGAGGCAGATTGATGGTATTGGATCAGGAGTCAATTTTCAGTTGAACAACCAACATAAATTCAACATTC
TAATCTTATATTCACTACCAGAAAGGAAAGAGATAGAGCAAGGGAAGAGCATAAAGTGCAGTTAACA
GATGTTGAGTCGACACAATGAAGGTGATGATCAACAAGGAAGCCGGTACAGTGTGATTCCACAGATTCAA
AAAGTGTGTGAAGTTGTAGATGGGTTTCATCTATGTTGCAAATGCTGAAGCTCATAAAAGTAAGTACTCAT
ATGTACATTTT

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC222474 representing NM_033484 Red=Cloning site Green=Tags(s) MAGSEPRSGTNSPPPPFSDWGRLEAAILSGWKTFWQSVSKERVARTTSREEVDEAASTLTRLPIDVQLYI LSFLSPHDLCLGSTNHYYNETVRDPIWRYFLLRDLPSSVDWKSLEILKKPISEVTDGAFFDYM AVYRMCCPYTRASKSSRPMYGAVTSFLHSLIIQNEPRFAMFGPGLEELNTSLVLSLMSSEELCPTAGLP QRQIDGIGSGVNFQLNNQHKFNILILYSTTRKERDRAREEHTSAVNKMF SRHNEGDDQQGSRYVIPQIQ KVCEVVDGFIYVANA EAHKSKYSYVHF SGPTRTRRLEQKLI SEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_033484
ORF Size:	921 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_033484.1, NP_277019.1</u>
RefSeq Size:	1746 bp
RefSeq ORF:	924 bp
Locus ID:	26272
UniProt ID:	<u>Q9UKT5</u>
Cytogenetics:	5p13.1
MW:	34.8 kDa