

## Product datasheet for **RC222503L3V**

### SLC25A24 (NM\_213651) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC25A24 (NM_213651) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC25A24
Synonyms:	APC1; SCAMC-1; SCAMC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_213651
ORF Size:	1374 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222503).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_213651.1</a>
RefSeq Size:	3344 bp
RefSeq ORF:	1377 bp
Locus ID:	29957
UniProt ID:	<a href="#">Q6NUK1</a>
Cytogenetics:	1p13.3
Protein Families:	Druggable Genome
MW:	51.2 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a carrier protein that transports ATP-Mg exchanging it for phosphate. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2012]