

## Product datasheet for **RC222535L4V**

### LEPRE1 (P3H1) (NM\_022356) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	LEPRE1 (P3H1) (NM_022356) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LEPRE1
Synonyms:	GROS1; LEPRE1; OI8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022356
ORF Size:	2208 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222535).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022356.2</a>
RefSeq Size:	2993 bp
RefSeq ORF:	2211 bp
Locus ID:	64175
UniProt ID:	<a href="#">Q32P28</a>
Cytogenetics:	1p34.2
Domains:	2OG-Fell_Oxy, P4Hc
Protein Families:	Secreted Protein



[View online »](#)

**MW:** 83.2 kDa

**Gene Summary:** This gene encodes an enzyme that is a member of the collagen prolyl hydroxylase family. These enzymes are localized to the endoplasmic reticulum and their activity is required for proper collagen synthesis and assembly. Mutations in this gene are associated with osteogenesis imperfecta type VIII. Three alternatively spliced transcript variants encoding different isoforms have been described. Other variants may exist, but their biological validity has not been determined. [provided by RefSeq, Aug 2011]