

## Product datasheet for **RC222544L3V**

### Neuroigin 2 (NLGN2) (NM\_020795) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Neuroigin 2 (NLGN2) (NM_020795) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NLGN2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020795
ORF Size:	2505 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222544).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020795.2</a>
RefSeq Size:	4642 bp
RefSeq ORF:	2508 bp
Locus ID:	57555
UniProt ID:	<a href="#">Q8NFZ4</a>
Cytogenetics:	17p13.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs)
MW:	90.6 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses. [provided by RefSeq, Jul 2008]