

Product datasheet for **RC222561L3V**

NFIX (NM_002501) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NFIX (NM_002501) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NFIX
Synonyms:	CTF; MRSHSS; NF-I/X; NF1-X; NF1A; SOTOS2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002501
ORF Size:	1323 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC222561).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002501.2
RefSeq Size:	5576 bp
RefSeq ORF:	1326 bp
Locus ID:	4784
UniProt ID:	Q14938
Cytogenetics:	19p13.13
Protein Families:	Transcription Factors
MW:	48.7 kDa


[View online »](#)

Gene Summary:

The protein encoded by this gene is a transcription factor that binds the palindromic sequence 5'-TTGGCNNNNNGCCAA-3 in viral and cellular promoters. The encoded protein can also stimulate adenovirus replication in vitro. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2012]